SMART **FITBODY** 2

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BRIEF INTRODUCTION

The product is a smart body fat scale with high-precision strain gauge sensors, which can help you store and keep track of your weight on a smartphone. This scale can transmit your weight data to your smartphone via wireless transmission technologies, as well as make graphs of your weight results, from which you can keep track of the trend anytime. By downloading the App on your phone and using it together with our scale, you will start your first step of weight control.

NOTE

This product can directly measure weight and perform body composition analysis on the personal information you input.

The body composition information provided by this product is for reference only and is not intended as medical advice. It is not recommended for use by children in their formative years, pregnant women, or people with artificial dialysis or heart bypass surgery, etc.

BATTERY INSTALLATION

It's powered by AAA batteries, please install as shown below:

First Use: Remove the insulation sheet before using (Note: if the battery is included).

Change battery: Open the battery cover, take out the used batteries, then install the new ones properly as directed shown below.

Note:

Please use the right battery. Remove battery if the scale leaves unused for more than 3 months.

The battery and related pictures are for reference only, please refer to the actual product.

HOW-TO-USE

Step 1: Place the scale properly and download the APP.

Place the scale on hard and flat floor.

- (1) Scan the QR code above to download and install "iMoni Health" APP.
or (2) Search for "iMoni Health" APP in "App Store" or "Google Play", then download and install it.

Step 2: Register and bind the device

Turn on the Bluetooth on your phone, open the APP, and follow the guidelines to register and bind the device:

The display is lighted up

Step gently to turn on the scale

Connect the scale with App

Warm reminder:

Personal body information is only used for body composition analysis, not for other purposes. To ensure the accuracy of the results, please enter accurate personal information.

Step 3: Measure and view the measurement data

Place the scale on a hard and flat floor. Stand on the scale in your bare feet. After the weight reading has stabilised and is flashing, the data transmission to the cloud via the mobile phone will begin. After a few seconds of successful transmission and display of weight data, the scale will automatically shut down. The APP displays and stores all index parameters and history.

Place the scale on hard and flat floor.

Keep your bare feet on the electrodes of scale.

Check the measurement data on App.

Tips:

- To make a proper measurement, please keep your bare feet on the electrodes of scale until the measuring is finished.
- If the APP is not connected for the first time of measuring, only weight can be measured, but no body fat. Please connect the APP before using.
- App is kept updating, so it may differ from showed on manual. For the final and accurate information, please check on App.

MULTIPLE USERS SUPPORTED

The scale can be connected with multiple users. Final user amount depends on the product.

DISPLAY ICONS

Low Battery Indication

The battery power is running low.

Calibration

Error exists, please re-test to get the correct result.

Over-load Indication

The weighing subject on the platform exceeds the maximum capability of scale. Please step off to avoid damage.

ABOUT THE UNIT

- It is suggested to set the unit on App.
- Or you can change the unit by the UNIT switch button on bottom of the scale(If the scale is with such button).
- The factory unit setting differs from different sales destination. Please refer to the actual product.

Back of the scale

ADVICE FOR USE & CARE

1. To ensure the accuracy of the data, it is recommended that each measurement be carried out at the same time of day. Don't use in strong electro-magnetic condition and field with intense shake .
2. When weighing after moving the scale or a long period without using it, please first press the scale surface to be zero and let it automatically shut down before turning on the scale.
3. Stand on the scale steadily and touch the electrodes on the platform with your bare feet. Tip over will lead to a fall when you step on the edge of one side of the scale. Keep standing firmly on scale when in use.
4. Do not use the scale in humid, overheat environment.(Temperature range from 5℃ to 35℃)
5. When the battery runs out, please remove the battery in time to prevent leakage of the battery lead scale body corrosion;
6. Do not drop, shock or strike the scale as it is a precise instrument.
7. Keep the scale clean. Do not use acid cleaning agent.
8. If the product does not work properly, please check that the batteries are inserted correctly or that the battery power is low. If necessary, reinstall the batteries and restart after 3 seconds or replace with new batteries.

FAQ

The software can't obtain weighing results from Smart scale.

- Please confirm the connection sign on the smart device to know the connection status.
- Confirm the Bluetooth is on in your device and the Smart scale is on within the valid range of signal coverage.
- Check whether the battery power is running low.
- Relaunch the application.
- Turn off the Bluetooth and on again, then bridge the Bluetooth connection between scale and smart phone again.

Fail to connect to your device

- Please make sure only one device is setting up connection with the Smart scale.
- Please make sure your device is equipped with Bluetooth 4.0 or above.

LCD has no response

- If it's powered by AAA batteries, please first check the batteries and replace them if necessary.
- Repairs may only be carried out by our authorized customer service or dealers.

ELDICO[™]
sport

In conformity with the
alternative packing
management system

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